








# EXTRAORAL | J



ITEM	DESCRIPTION	CATALOG NO.																																												
	<p><b>Extraoral Facebow (Pro-Bow™)</b> This facebow is comfortable for the patient, plus it is easy to adjust and durable. The bow is made of solid wire with no hollow wire inserted. Hygienic, corrosion resistant, and more visually attractive. Wire diameter: .045; Angulation: 0°. Vertical loops on inner bow are adjustable.</p> <table border="1" data-bbox="561 531 1235 810"> <thead> <tr> <th>INNER ARCH</th> <th>SIZE NUMBER</th> <th>MM LENGTH</th> <th>CATALOG NO.</th> </tr> </thead> <tbody> <tr> <td>.045</td> <td>1</td> <td>83</td> <td>425-011A</td> </tr> <tr> <td>.045</td> <td>2</td> <td>90</td> <td>425-012A</td> </tr> <tr> <td>.045</td> <td>3</td> <td>97</td> <td>425-013A</td> </tr> <tr> <td>.045</td> <td>4</td> <td>104</td> <td>425-014A</td> </tr> <tr> <td>.045</td> <td>5</td> <td>111</td> <td>425-015A</td> </tr> <tr> <td>.045</td> <td>6</td> <td>Straight</td> <td>425-016A</td> </tr> </tbody> </table>	INNER ARCH	SIZE NUMBER	MM LENGTH	CATALOG NO.	.045	1	83	425-011A	.045	2	90	425-012A	.045	3	97	425-013A	.045	4	104	425-014A	.045	5	111	425-015A	.045	6	Straight	425-016A																	
INNER ARCH	SIZE NUMBER	MM LENGTH	CATALOG NO.																																											
.045	1	83	425-011A																																											
.045	2	90	425-012A																																											
.045	3	97	425-013A																																											
.045	4	104	425-014A																																											
.045	5	111	425-015A																																											
.045	6	Straight	425-016A																																											
	<p><b>Safety Neck Pads</b> These modules greatly reduce the risk of slingshot rebound from headgear. Reassembly after release is easy. The stainless steel coil spring allows the doctor to apply specific forces depending on in which hole he or she inserts the facebow.</p> <p>Neck Pad &amp; Safety Release with Stainless Steel Spring</p> <ul style="list-style-type: none"> <li>• SS Black Safety Neck Pad (Black Clip)</li> <li>• Black Pad Component (5 per pack)</li> </ul>	<p><b>426-002</b> <b>426-000</b></p>																																												
	<p><b>Velpro™ Safety Neck Pad</b> Soft <i>Velpro</i> Material offers maximum patient comfort, while greatly reducing slingshot rebound from headgear. Reassembly after release is easy. The stainless steel coil spring allows the doctor to apply specific forces depending on which hole he or she inserts the facebow into. Comes as a complete assembly only with black clip.</p> <p>Black Purple Pink Blue Teal Black Clip Component (10 per pack)</p> <p><b>Force of Clip Component</b></p> <table border="1" data-bbox="448 1446 1073 1818"> <thead> <tr> <th>INNER ARCH</th> <th>SIZE NUMBER</th> <th>MM LENGTH</th> <th>CATALOG NO.</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>291</td> <td>10</td> <td>0.65</td> </tr> <tr> <td>2</td> <td>414</td> <td>15</td> <td>0.98</td> </tr> <tr> <td>3</td> <td>524</td> <td>19</td> <td>1.17</td> </tr> <tr> <td>4</td> <td>613</td> <td>22</td> <td>1.37</td> </tr> <tr> <td>5</td> <td>730</td> <td>26</td> <td>1.63</td> </tr> <tr> <td>6</td> <td>875</td> <td>31</td> <td>1.95</td> </tr> <tr> <td>7</td> <td>966</td> <td>35</td> <td>2.16</td> </tr> <tr> <td>8</td> <td>1104</td> <td>39</td> <td>2.47</td> </tr> <tr> <td>9</td> <td>1196</td> <td>43</td> <td>2.67</td> </tr> <tr> <td>10</td> <td>1342</td> <td>49</td> <td>3.04</td> </tr> </tbody> </table>	INNER ARCH	SIZE NUMBER	MM LENGTH	CATALOG NO.	1	291	10	0.65	2	414	15	0.98	3	524	19	1.17	4	613	22	1.37	5	730	26	1.63	6	875	31	1.95	7	966	35	2.16	8	1104	39	2.47	9	1196	43	2.67	10	1342	49	3.04	<p><b>426-070</b> <b>426-071</b> <b>426-072</b> <b>426-073</b> <b>426-074</b> <b>426-020</b></p>
INNER ARCH	SIZE NUMBER	MM LENGTH	CATALOG NO.																																											
1	291	10	0.65																																											
2	414	15	0.98																																											
3	524	19	1.17																																											
4	613	22	1.37																																											
5	730	26	1.63																																											
6	875	31	1.95																																											
7	966	35	2.16																																											
8	1104	39	2.47																																											
9	1196	43	2.67																																											
10	1342	49	3.04																																											

ITEM	DESCRIPTION	CATALOG NO.
	<p><b>Face Mask (Reverse Pull)</b> It is multi-adjustable and features a forehead and chin rest with dual air vents for maximum patient comfort and air circulation. Chin cup can be adjusted for a fixed position or for full range of motion. It is designed to distribute pressure evenly, allowing for significant range of motion without side-to-side rotation. The few bonus features are:</p> <ul style="list-style-type: none"> <li>• Maximum patient comfort</li> <li>• High impact molded plastic construction</li> <li>• Lightweight and fully adjustable</li> <li>• Perfect for patients that wear glasses</li> <li>• Includes two sets of replacement pads and one hexagon key</li> </ul> <p>Blue Red Teal Pearl</p>	<p>426-055 426-056 426-057 426-058</p>
	<p><b>Extraoral Elastics</b> Available only in amber, 1/2" thickness (1000 per pack) 8 oz 16 oz</p>	<p>408-073 408-074</p>
	<p><b>Safety High-Pull Head Caps</b> Extra heavy duty construction for sturdy dependability. Easy to adjust.</p> <p><b>SS Black High Pull Cap</b></p> <ul style="list-style-type: none"> <li>• Black Clip</li> </ul> <p><b>Safety High-Pull Head Cap Components</b></p> <ul style="list-style-type: none"> <li>• Black Caps (5 per pack)</li> <li>• Black Clips (10 per pack)</li> </ul>	<p>426-012 426-010 426-020</p>